

Human DAA Recombinant Protein

Catalog Number: PROTP14920

Overview

| | |
|--------------|--|
| Product Name | Human DAA Recombinant Protein |
| Description | Human DAA Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-347aa. Application(s): SDS-PAGE. |
| Size | Starting from 10ug |
| Tag | His-Tag |
| Form | Liquid |
| Source | E. coli |
| Formulation | 20mM Tris-HCl buffer (pH 8.0) containing 20% glycerol 1mM DTT |

Concentration

0.5mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Human, P14920, 1-347aa; MRVVVIGAGV IGLSTALCIH ERYHSVLQPL DIKVYADRFT PLTTTDVAAG LWQPYLSDPN NPQEADWSQQ
TFDYLLSHVH SPNAENLGLF LISGYNLFHE AIPDPSWKDT VLGFRKLTPR ELDMPDYGY GWFHTSLILE GKNYLQWLTE RLTERGVKFF
QRKVESFEEV AREGADVIVN CTGWWAGALQ RDPLLQPGRG QIMKVDAPWM KHFILTHDPE RGIYNSPYII PGTQVTLLGG IFQLGNWSEL
NNIQDHNTIW EGCCRLEPTL KNARIIGERT GFRPVRPQIR LEREQLRTGP SNTVEIHNYG HGGYGLTIHW GCALEAAKLF GRILEEKKLS RMPPSHL

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human DAA Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.