

## Human ETS2 Recombinant Protein

Catalog Number: PROTP15036

### Overview

Product Name	Human ETS2 Recombinant Protein
Description	Human ETS2 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-469aa. Application(s): SDS-PAGE, Denatured. Denatured protein is most likely not the best option for functional studies. It is better suited for Western Blot (WB) or imaging assays.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 10% glycerol, 0.4M urea

### Concentration

1mg/ml (determined by Bradford assay)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 85% by SDS-PAGE

### Amino Acid Sequence

Human, P15036, 1-469aa; MNDFGIKNMD QVAPVANSYR GTLKRQPAFD TFDGSLFAVF PSLNEEQTLQ EVPTGLDSIS HDSANCELPL  
LTPCSKAVMS QALKATFSGF KKEQRRLGIP KNPWLWSEQQ VCQWLLWATN EFSLVNVLNQ RFGMNGQMLC NLGKERFLEL APDFVGDILW  
EHLEQMIKEN QEKTEDQYEE NSHLTSVPHW INSNTLGFGE EQAPYGMQTQ NYPKGGLLDS MCPASTPSVL SSEQEFQMFP KSRLSSVSVT  
YCSVSQDFPG SNLNLTTNNS GTPKDHDSP E NGADSFESSD SLLQSWNSQS SLLDVQRVPS FEFEDDCSQ SLCLNKPTMS FKDYIQERSD  
PVEQKGKPVIP AAVLAGFTGS GPIQLWQFLL ELLSDKSCQS FISWTGDGWE FKLADPDEVA RRWGKRKNKP KMNYEKLSRG LRYYYDKNII  
HKTS GKRYVY RFVCDLQNLGFTPEELHAI LGVQPDTE

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human ETS2 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.