

## Human recombinant CD19 protein

Catalog Number: PROTP15391-1

### Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Overview

|              |                                                                                                                                                                                                  |
|--------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name | Human recombinant CD19 protein                                                                                                                                                                   |
| Description  | Human recombinant CD19 protein expressed in Baculovirus with hlgG-His-Tag. Sequence domain: 21-291aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method). |
| Size         | Starting from 10ug                                                                                                                                                                               |
| Tag          | hlgG-His-Tag                                                                                                                                                                                     |
| Form         | Liquid                                                                                                                                                                                           |
| Source       | Baculovirus                                                                                                                                                                                      |
| Formulation  | Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol                                                                                                                                       |

### Concentration

0.25mg/ml (determined by absorbance at 280nm)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 85% by SDS-PAGE

### Amino Acid Sequence

Human, P15391, 21-291aa; EEPLVVKVEE GDNAVLQCLK GTSDGPTQQL TWSRESPLKP FLKLSLGLPG LGIHMRLAI WLFIFNVSQQ  
MGGFYLCQPG PPSEKAWQPG WTVNVEGSGE LFRWNVSDLG GLGCGLKNRS SEGSSPSGK LMSPKLYVWA KDRPEIWEGE PPCLPPRDSL  
NQSLSQDLTM APGSTLWLSC GVPPDSVSRG PLSWTHVHPK GPKSLLSLEL KDDRPARDMW VMETGLLLPR ATAQDAGKYY CHRGNLTMSF  
HLEITARPVL WHWLLRTGGW K

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

### [Human recombinant CD19 protein](#)

For Research Use Only. Not for use in diagnostic procedures.