

Human 15-PGDH/HPGD Recombinant Protein

Catalog Number: PROTP15428

Overview

Product Name	Human 15-PGDH/HPGD Recombinant Protein
Description	Human 15-PGDH/HPGD Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-266aa. Application(s): SDS-PAGE.
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 20% glycerol, 0.1M NaCl, 1mM DTT

Concentration

0.5mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Human, P15428, 1-266aa; MHVNGKVALV TGAAQGIGRA FAEALLLKGA KVALVDWNLE AGVQCKAALD EQFEPQKTLF IQCDVADQQQ
LRDTRFRKVV D HFGRLDILVN NAGVNNEKNW EKTQLINLVS VISGTYLGLD YMSKQNGGEG GIINMSSLA GLMPVAQQPV YCASKHGIVG
FTRSAALAN LMNSGVRLNA ICPGFVNTAI LESIEKEENM GQYIEYKDHI KDMIKYYGIL DPPLIANGLI TLIEDDALNG AIMKITTSKG IHFQDYDTTP
FQAKTQ

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human 15-PGDH/HPGD Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.