

## Human Vav-1 Recombinant Protein

Catalog Number: PROTP15498

### Overview

Product Name	Human Vav-1 Recombinant Protein
Description	Human Vav-1 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 189-565aa. Application(s): SDS-PAGE.
Size	Starting from 4ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 0.2M NaCl, 30% glycerol, 2mM DTT

### Concentration

0.25mg/ml (determined by Bradford assay)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 85% by SDS-PAGE

### Amino Acid Sequence

Human, P15498, 189-565aa; MGSSHHHHHH SGLVPRGSH MGSMT EYDKR CCCLREIQQT EEKYDTLGS IQQHFLKPLQ RFLKPQDIEI  
IFINIEDLLR VHTHFLKEMK EALGTPGAAN LYQVFIKYKE RFLVYGRYCS QVESASKHLD RVAAAREDVQ MKLEEC SQRA NNGRFTLRDL  
LMVPMQRVLK YHLLLQELVK HTQEAMEKEN LRLALDAMRD LAQCVNEVKR DNETLRQITN FQLSIENLDQ SLAHYGRPKI DGELKITSVE  
RRSKMDRYAF LLDKALLICK RRGDSYDLKD FVNLHSFQVR DDSSGDRDNK KWSHMFLLE DQGAQGYELF FKTRCLKKKW MEQFEMAISN  
IYPENATANG HDFQMFSFEE TTSCACQML LRGTIFYQGYR CHR CRASAHK ECLGRVPPCG

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Human Vav-1 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.