

Human recombinant CD1d protein

Catalog Number: PROTP15813-1

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Human recombinant CD1d protein
Description	Human recombinant CD1d protein expressed in Baculovirus with His-Tag. Sequence domain: 20-301aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

0.25mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 85% by SDS-PAGE

Amino Acid Sequence

Human, P15813, 20-301aa; EVPQRLFPLR CLQISSFANS SWTRTDGLAW LGELQTHSWS NDSDTVRSRK PWSQGTFSQD QWETLQHIFR
VYRSSFTRDV KEFAKMLRLS YPLELQVSAG CEVHPGNASN NFFHVAFQ GK DILSFQGT SW EPTQEAPLWV NLAIQVLNQD KWTRETQVWL
LNGTCPQFVS GLLES GKSEL KKQVKPKAWL SRGSPGPGR LLLVCHVSGF YPKPVWVKWM RGEQEQQGTQ PGDILPNADE TWYLRATL DV
VAGEAAGLSC RVKHSSLEGQ DIVLYWGGSY TS

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

[Human recombinant CD1d protein](#)

For Research Use Only. Not for use in diagnostic procedures.