

Mouse HGF R/c-MET Recombinant Protein

Catalog Number: PROTP16056

Overview

Product Name	Mouse HGF R/c-MET Recombinant Protein
Description	Mouse HGF R/c-MET Recombinant Protein expressed in Baculovirus with hIgG-His-tag. Sequence domain: 25-931aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	hIgG-His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

0.25mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 85% by SDS-PAGE

Amino Acid Sequence

Mouse, P16056, 25-931aa; ECKEALVKSE MNVNMKYQLP NFTAETPIQN VVLHGHHIYL GATNYIYVLN DKDLQKVSEF KTGPVLEHPD
CLPCRDCSSK ANSSGGVWKD NINMALLVDT YYDDQLISCG SVNRGTCQRH VLPPDNSADI QSEVHCMFSP EESGQCPCD VVSALGAKVL
LSEKDRFINF FVGNTINSSY PPGYSLHSIS VRRLKETQDG FKFLTDQSYI DVLPEFQDSY PIKYIHAFES NHFIYFLTVQ KETLDAQTFH TRIIRFCSVD
SGLHSYMEMP LECILTEKRR KRSTREEVFN ILQAAYVSKP GANLAKQIGA SPSDDILFGV FAQSKPDSAE PVNRSVCAAF PIKYVNDFFN
KIVNKNNVRC LQHFYGNHE HCFNRTLLRN SSGCEARSDE YRTEFTTALQ RVDLFMGRLN QVLLTSISTF IKGDLTIANL GTSEGRFMQV
VLSRTAHLTP HVNFLDSDP VSPEVIVEHP SNQNGYTLVV TGKKITKIPL NGLGCGHFQS CSQCLSAPYF IQCGWCHNQC VRFDECPST
WTQEICLPAV YKVFPTSAPL EGGTVLTICG WDFGFRKNNK FDLRKTQVLL GNESTLTLT ESTTNTLKCT VGPAMSEHFN VSVIISNSRE
TTQYSAFSYV DPVITSISPR YGPQAGGTLT TLTKGKYLNSG NSRHISIGGK TCTLKSVSDS ILECYTPAQT TSDEFVVKLK IDLANRETSS FSYREDPVVY
EIHPTKSFIS GGSTITGIGK TLNSVSLPKL VIDVHEVGVN YTVACQHRSN SEIICCTTPS LKQLGLQLPL KTKAFFLLDG ILSKHFDLTY VHNPFVPEFE
KPVMIISIGNE NVVEIKGNNI DPEAVKGEVL KVGNGQSCESL HWHSGAVLCT VPSDLLKLN ELNIEWKQAV SSTVLGKIV QPDQNF

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Mouse HGF R/c-MET Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.