

Human CD43 Recombinant Protein

Catalog Number: PROTP16150

Overview

Product Name	Human CD43 Recombinant Protein
Description	Human CD43 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 20-253aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 10% glycerol

Concentration

1mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 85% by SDS-PAGE

Amino Acid Sequence

Human, P16150, 20-253aa; STTAVQT PTSGEPLVST SEPLSSKMYT TSITSDPKAD STGDQTSALP PSTSINEGSP LWTSIGASTG SPLPEPTYQ
EVSIKMSSVP QETPHATSHP AVPITANSLG SHTVTGGTIT TNSPETSRT SGAPVTTAAS SLETSRGTSG PPLTMATVSL ETSKGTSGPP VTMATDSLET
STGTTGPPVT MTTGSLEPSS GASGPQVSSV KLSTMMSPTT STNASTVPFR NPDENSR

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human CD43 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.