

Human PDGF R alpha Recombinant Protein

Catalog Number: PROTP16234

Overview

Product Name	Human PDGF R alpha Recombinant Protein
Description	Human PDGF R alpha Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 24-524aa. Application(s): SDS-PAGE.
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 20% glycerol, 0.1M NaCl, 1mM DTT

Concentration

0.5mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 85% by SDS-PAGE

Amino Acid Sequence

Human, P16234, 24-524aa; MGSSHHHHHH SGLVPRGSH MQLSLPSILP NENEKVVQLN SSFSLRCFGE SEVSWQYPMS EEESDVEIR
NEENNSGLFV TVLEVSSASA AHTGLYTCYY NHTQTEENEL EGRHIYIYVP DPDVAFVPLG MTDYLVIVED DDSAIIPCRT TDPETPVTLH
NSEGVVPASY DSRQGFNGTF TVGPYICEAT VKGKKFQTIP FNVYALKATS ELDLEMEALK TVYKSGETIV VTCAVFNNEV VDLQWTPGE
VKGKGITMLE EIKVPSIKLV YTLTVPEATV KDSGDYECOA RQATREVKEM KKVTISVHEK GFIEIKPTFS QLEAVNLHEV KHFVVEVRAY
PPPRISWLKN NLTLIENLTE ITTDVEKIQE IRYRSKLLI RAKEEDSGHY TIVAQNEDAV KSYTFELLTQ VPSSILDLDV DHHGSTGGQT VRCTAEGTPL
PDIEWMICKD IKKCNNETSW TILANNVSN I TEIHSRDRS TVEGRVTFK VEETIAVRCL AKNLLGAENR ELKLVAPTLR SE

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human PDGF R alpha Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.