

Human CD31/PECAM-1 Recombinant Protein

Catalog Number: PROTP16284

Overview

Product Name	Human CD31/PECAM-1 Recombinant Protein
Description	Human CD31/PECAM-1 Recombinant Protein expressed in Baculovirus with His-tag. Sequence domain: 28-601aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method). Bioactivity: Measured by the ability of the immobilized protein to support the adhesion of Jurkat human acute T cell leukemia cells. When cells are added to human CD31 coated plates 5 ug/ml. This effect is more to 40%.
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

0.25mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Human, P16284, 28-601aa; QENSFTINSV DMKSLPDWTV QNGKNLTLQC FADVSTTSHV KPQHQLMLFYK DDVLFYNISS MKSTESYFIP
EVRIYDSGT Y KCTVIVNNKE KTTAEYQVLV EGVSPRVTL DKKEAIQGGI VRVNCVPEE KAPIHFTIEK LELNEKMKV KL KREKNSRDQN
FVILEFPVEE QDRVLSFRCQ ARIISGIHMV TSESTKSELV TVTESFSTPK FHISPTGMIM EGAQLHIKCT IQVTHLAQEF PEIIIQKDKA IVAHNRHGK
AVYSVMAMVE HSGNYTCKVE SSRISKVSSI VVNITELFSK PELESSFTHL DQGERLNLSC SIPGAPPANF TIQKEDTIVS QTQDFTKIAS KSDSGTYICT
AGIDKVVKKS NTVQIVVCEM LSQPRISYDA QFEVIKGGTI EVRCESISGT LPISYQLLKT SKVLENSTKN SNDAVFKDN PTEDVEYQCV
ADNCHSHAKM LSEVLRVKVI APVDEVQISI LSSKVVESGE DIVLQCAVNE GSGPITYKFY REKEGKPFYQ MTSNATQAFW TKQKASKEQE
GEYYCTAFNR ANHASSVPRS KILTVRVILA PWKKVEHHHH HH

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human CD31/PECAM-1 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.