

Rat Growth Hormone R/GHR Recombinant Protein

Catalog Number: PROTP16310

Overview

Product Name	Rat Growth Hormone R/GHR Recombinant Protein
Description	Rat Growth Hormone R/GHR Recombinant Protein expressed in Baculovirus with His-tag. Sequence domain: 19-265aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

0.5mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles.

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Rat, P16310, 19-265aa; FPGSGATPAT LGKASPVLR INPSLRESSS GKPRFTKCRS PELETFCSCYW TEGDDHNLKV PGSIQLYYAR RIAHEWTPPEW
KECPDYVSAG ANSCYFNSSY TSIWIPYCIK LTTNGDLLKE KCFTVDEIVQ PDPIGLNWT LLNISLPGIR GDIQVSWQPP PSADVLKGWI ILEYEIQYKE
VNETKWKTMS PIWSTSVPLY SLRLDKEHEV RVRSRQRSFE KYSEFSEVLR VTFPQMDTLA ACEEDFRLEH HHHHH ?/description_after_attributes?

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Rat Growth Hormone R/GHR Recombinant Protein