

CTLA4 (NM_005214) Human Recombinant Protein

Catalog Number: PROTP16410

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	CTLA4 (NM_005214) Human Recombinant Protein
Description	Purified recombinant protein of Human cytotoxic T-lymphocyte-associated protein 4 (CTLA4), transcript variant 1, with C-terminal human Fc and his tags, expressed in human cells
Size	20 µg
Tag	C-Fc/His
Form	Frozen Solution in PBS Buffer
Source	HEK293T
Formulation	PBS, pH7.4, 10% glycerol

Concentration

>50 ug/mL as determined by microplate BCA method

Storage

Store at -80°C. Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles. (Ship on dry ice.)

Purity

> 80% as determined by SDS-PAGE and Coomassie blue staining

Amino Acid Sequence

A DNA sequence encoding the extracellular domain (Ala37-Phe162) of human CTLA4

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

CTLA4 (NM_005214) Human Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.