

Human recombinant CD36/SR-B3 protein

Catalog Number: PROTP16671-1

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Human recombinant CD36/SR-B3 protein
Description	Human recombinant CD36/SR-B3 protein expressed in Baculovirus with hlgG-His-Tag. Sequence domain: 30-439aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	hlgG-His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

0.5mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 85% by SDS-PAGE

Amino Acid Sequence

Human, P16671, 30-439aa; ADLGDLIIQK TIKKQVVLEE GTIAFKNWVK TGTEVYRQFW IFDVQNPQEV MMNSSNIQVK QRGPTYRVR
FLAKENVTDQ AEDNTVSFLQ PNGAIFEPSL SVGTEADNFT VLNLAVAAAS HIYQNQFVQM ILNSLINKSK SSMFQVRTLR ELLWGYRDPF
LSLVPPVTT TVGLFYPYNN TADGVYKVFN GKDNISKVAI IDTYKGKRNLSY WESHCDMI NGTDAASFPP FVEKSQVLQF FSSDICRSIY
AVFESDVNLK GIPVYRFVLP SKAFASPVEN PDNYCFCTEK IISKNCTSYG VLDISKCKEG RPVYISLPHF LYASPDVSEP IDGLNPNEEE HRTYLDIEPI

TGFTLQFAKR LQVNLLVKPS EKIQLKLNK RNYIVPILWL NETGTIGDEK ANMFRSQVTG KINLEPKSCD KTHTCPPCPA PELLGGPSVF
LFPPKPKDTL MISRTPEVTC VVVDVSHEDP EVKFNWYVDG VEVHNAKTKP REEQYNSTYR VVSVLTVLHQ DWLNGKEYKC KVSNAKALPAP
IEKTISKAKG QPREPQVYTL PPSRDELTKN QVSLTCLVKG FYPSDIAVEW ESNQGPENNY KTTTPVLDSG GSFFLYSKLT VDKSRWQQGN
VFSCSVMHEA LHNHYTQKSL SLSPGKHHHH HH

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

[Human recombinant CD36/SR-B3 protein](#)

For Research Use Only. Not for use in diagnostic procedures.