

Mouse Lamp2 Recombinant Protein

Catalog Number: PROTP17047

Overview

| | |
|--------------|--|
| Product Name | Mouse Lamp2 Recombinant Protein |
| Description | Mouse Lamp2 Recombinant Protein expressed in Baculovirus with His-tag. Sequence domain: 26-379aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method). |
| Size | Starting from 10ug |
| Tag | His-Tag |
| Form | Liquid |
| Source | Baculovirus |
| Formulation | Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol |

Concentration

0.5mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Mouse, P17047, 26-379aa; LIVNLTDSKG TCLYAEWEMN FTITYETTQ TNKTITIAVP DKATHDGSSC GDDRNSAKIM IQFGFAVSWA
VNFTKEASHY SIHDIVLSYN TSDSTVPGA VAKGVHTVKN PENFKVPLDV IFKCNVLTLY NLTPVVQKYW GIHLQAFVQN GTVSKNEQVC
EEDQTPTTVA PIIHTTAPST TTTLTPTSTP TPTPTPTPTV GNYSIRNGNT TCLLATMGLQ LNITEEKVPF IFNINPATTN FTGSCQPQSA QLRNNSQIK
YLDIFFAVKV EKR FYLKEVN VYMYLANGSA FNISNKNLSF WDAPLGSSYM CNKEQVLSVS RAFQINTFNL KVQPFNVTKG QYSTAQDCSA DEDN

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Mouse Lamp2 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.