

Human IFN-alpha/beta R1 Recombinant Protein

Catalog Number: PROTP17181

Overview

Product Name	Human IFN-alpha/beta R1 Recombinant Protein
Description	Human IFN-alpha/beta R1 Recombinant Protein expressed in Baculovirus with His-tag. Sequence domain: 28-436aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 20% glycerol

Concentration

0.5mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Human, P17181, 28-436aa; KNLKSPQKVE VDIIDDNFIL RWNRSDES VG NVTFSFDYQK TGMDNWIKLS GCQNITSTKC NFSSLKLN VY EEIKLRIRAE KENTSSWYEV DSFTPF RKAQ IGPPEVHLEA EDKAI VIHIS PGT KD SVMWA LDGLSFTYSL VIWKNSSGVE ERIENIYSRH KIYKLS PETT YCLKVKAALL TSWKIGVYSP VHCIKTTVEN ELPPPENIEV SVQNQNYVLK WDYTYANMTF QVQWLHAF LK RNP GNHLYKW KQIPDCENVK TTQCVFPQNV FQKGIYLLRV QASDGNNTSF WSEEIKFDTE IQAFL LPPVF NIRSLSDSFH IYIGAPKQSG NTPVIQDYPL IYEIIFWENT SNAERKIEK KTDVTVPNLK PLTVYCVKAR AHTMDEKLNK SSVFSDAVCE KTKPGNTSK

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human IFN-alpha/beta R1 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.