

Human PKA C alpha Recombinant Protein

Catalog Number: PROTP17612

Overview

Product Name	Human PKA C alpha Recombinant Protein
Description	Human PKA C alpha Recombinant Protein expressed in Baculovirus with GST-tag. Sequence domain: 1-351aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	GST-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

0.25mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 85% by SDS-PAGE

Amino Acid Sequence

Human, P17612, 1-351aa; MGN AAAAKKGSEQ ESVKEFLAKA KEDFLKKWES PAQNTAHL DQ FERIKTLGTG SFGRVMLVKH KETGNHYAMK
ILDKQKVVKL KQIEHTLNEK RILQAVNFPF LVKLEFSFKD NSNLYMMEY VPGGEMFSL RRIGRFSEPH ARFYAAQIVL TFEYLHSLDL
IYRDLKPENL LIDQQGYIQV TDFGFAKRVK GRTWTLCGTP EYLAPEIILS KGYNKAVDWW ALGVLIYEMA AGYPPFFADQ PIQIYEKIVS
GKVRFP SHFS SDLKDLLRNL LQVDLTKRFG NLKNGVNDIK NHHKWFATTDW IAIYQRKVEA PFIKFKGPG DTSNFD DYE E EIRVSINEK CGKEFSEF

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human PKA C alpha Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.