

Human HLA-G Recombinant Protein

Catalog Number: PROTP17693

Overview

Product Name	Human HLA-G Recombinant Protein
Description	Human HLA-G Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 25-308aa. Application(s): SDS-PAGE.
Size	Starting from 4ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 20% glycerol

Concentration

0.25mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, P17693, 25-308aa; GSHSM RYFSAAVSRP GRGEPRIAM GYVDDTQFVR FDSDSACPRM EPRAPWVEQE GPEYWEEETR NTKAHAQTDR
MNLQTLRGYY NQSEASSHTL QWMIGCDLGS DGRLLRGYEQ YAYDGKDYLA LNEDLRSWTA ADTAAQISKR KCEAANVAEQ RRAYLEGTCV
EWLHRYLENG KEMLQRADPP KTHVTHHPVF DYEATLRCWA LGFYPAEIL TWQRDGEDQT QDVELVETRP AGDGTQKWA AVVVPSGEEQ
RYTCHVQHEG LPEPLMLRWK QSSLPTIPI

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human HLA-G Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.