

e.coli mutY Recombinant Protein

Catalog Number: PROTP17802

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	e.coli mutY Recombinant Protein
Description	E. coli mutY Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-350aa. Application(s): SDS-PAGE.
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 0.1M NaCl, 20% glycerol

Concentration

0.25mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

E.coli, P17802, 1-350aa; MGSSHHHHHH SGLVPRGSH MGSMQASQFS AQVLDWYDKY GRKTLPWQID KTPYKVLSE VMLQQTQVAT
VIPYFERFMA RFPTVTDLAN APLDEVHLHW TGLGYARAR NLHKAAQQA TLHGGKFPET FEEVAALPGV GRSTAGAILS LSLGKHFPIL
DGNVKRVLAR CYAVSGWPGK KEVENKLWSL SEQVTPAVGV ERFNQAMMDL GAMICTRSKP KCSLCPLQNG CAAAANNSWA LYPGKKPKQT
LPERTGYFLL LQHEDEVLLA QRPPSGLWGG LYCFPQFADE ESLRQWLAQR QIAADNLTQL TAFRHTFSHF HLDIVPMWLP VSSFTGCMDE
GNALWYNLAQ PPSVGLAAPV ERLQLQRTG APV

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

[e.coli mutY Recombinant Protein](#)

For Research Use Only. Not for use in diagnostic procedures.