

## Human recombinant LBP protein

Catalog Number: PROTP18428-1

### Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Overview

|              |                                                                                                                                                                                       |
|--------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name | Human recombinant LBP protein                                                                                                                                                         |
| Description  | Human recombinant LBP protein expressed in HEK293 with His-Tag. Sequence domain: 26-481aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method). |
| Size         | Starting from 10ug                                                                                                                                                                    |
| Tag          | His-Tag                                                                                                                                                                               |
| Form         | Liquid                                                                                                                                                                                |
| Source       | HEK293                                                                                                                                                                                |
| Formulation  | Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol                                                                                                                            |

### Concentration

0.25mg/ml (determined by absorbance at 280nm)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 95% by SDS-PAGE

### Amino Acid Sequence

Human, P18428, 26-481aa; ANPGLVARIT DKGLQYAAQE GLLALQSELL RITLPDFTGD LRIPHVGRGR YEFHSLNIHS CELLHSALRP VPGQGLSLSI  
SDSSIRVQGR WKVRKSFFKL QGSFDVSVKG ISISVNLLG SESSGRPTVT ASSCSSDIAD VEVDMSGDLG WLLNLFHNQI ESKFQKVLES  
RICEMIQKSV SSDLQPYLQT LPVTTEIDSF ADIDYSLVEA PRATAQMLEV MFKGEIFHRN HRSPVTLLAA VMSLP EEHNK M VYFAISDYV  
FNTASLVYHE EGYLNFSITD DMIPPDNIR LTTKSRPFV PRLARLYPNM NLELQGSVPS APLLNFSPGN LSVDPYMEID AFVLLPSSSK EPVFRLSVAT  
NVSATLTFNT SKITGFLKPG KVKVELKESK VGLFNAELLE ALLNYYILNT FYPKFNDKLA EGFLPLLKR VQLYDLGLQI HKDFLFLGAN VQYMRV

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

### Human recombinant LBP protein

For Research Use Only. Not for use in diagnostic procedures.