

Human recombinant IL-1Ra/IL1RN protein

Catalog Number: PROTP18510-4

Overview

Product Name	Human recombinant IL-1Ra/IL1RN protein
Description	Human recombinant IL-1Ra/IL1RN protein expressed in E.coli with Non-Tagged. Sequence domain: 26-177aa. Application(s): SDS-PAGE, Bioactivity. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	Non-Tagged
Form	Liquid
Source	E.coli
Formulation	Phosphate-Buffered Saline (pH7.4)

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Human, P18510, 26-177aa; MRPSGRKSSK MQAFRIWDVN QKTFYLRNNQ LVAGYLQGPV VNLEEKIDVV PIEPHALFLG IHGGKMCLSC VKSGDETRLQ LEAVNITDLS ENRKQDKRFA FIRSDSGPTT SFESAACPGW FLCTAMEADQ PVSLTNMPDE GVMVTKFYFQ EDE

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human recombinant IL-1Ra/IL1RN protein

For Research Use Only. Not for use in diagnostic procedures.