

Human RDBP Recombinant Protein

Catalog Number: PROTP18615

Overview

Product Name	Human RDBP Recombinant Protein
Description	Human RDBP Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-380aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 50% glycerol, 1mM DTT, 0.1mM PMSF

Concentration

0.25mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 85% by SDS-PAGE

Amino Acid Sequence

Human, P18615, 1-380aa; MLVIPPGLSE EEEALQKKFN KLKKKKKALL ALKKQSSST TSQGGVKRSL SEQPVMDTAT ATEQAKQLVK SGAISAIKAE
TKNSGFKRSR TLEGKLDPE KGPVPTFQPF QRSISADDDL QESSRRPQRK SLYESFVSSS DRLRELPGDG EEAEGPGAGD GPPRSFDWGY
EERSGAHSSA SPPRSRSRDR SHERNRDRDR DRERDRDRDR DRDRERDRDR DRDRDRDRER DRDRERDRDR DREGPFRRSD SFPERRAPRK
GNTLYVYGED MTPTLLRGAF SPFGNIIDLS MDPPRNCAFV TYEKMESADQ AVAELNGTQV ESVQLKVNIA RKQPMLDAAT GKSVWGLAV
QNSPKGCHRD KRTQIVYSDD VYKENLVDGF

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human RDBP Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.