

Human Syndecan-1/CD138 Recombinant Protein

Catalog Number: PROTP18827

Overview

Product Name	Human Syndecan-1/CD138 Recombinant Protein
Description	Human Syndecan-1/CD138 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 18-254aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 0.1M NaCl, 10% glycerol

Concentration

1mg/ml (determined by BCA assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, P18827, 18-254aa; QPALP QIVATNLPPE DQDGS GDDSD NFSGSGAGAL QDITLSQQTP STWKDTQLLT AIPTSPEPTG LEATAASTST LPAGEGPK EG EAVVLPEVEP GLTAREQEAT PRPRETTQLP TTHQASTTTA TTAQEPATSH PHRDMQPGHH ETSTPAGPSQ ADLHTPHTED GGPSATERAA EDGASSQLPA AEGSGEQDFT FETSGENTAV VAVEPDRRNQ SPVDQGATGA SQGLLDRKEV LG

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human Syndecan-1/CD138 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.