

Human Hexokinase 1 Recombinant Protein

Catalog Number: PROTP19367

Overview

Product Name	Human Hexokinase 1 Recombinant Protein
Description	Human Hexokinase 1 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-917aa. Application(s): SDS-PAGE, Enzyme Activity.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 10% glycerol

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 85% by SDS-PAGE

Amino Acid Sequence

Human, P19367, 1-917aa; MIAAQLLAYY FTELKDDQVK KIDKLYAMR LSETLIDIM TRFRKEMKNG LSRDFNPTAT VKMLPTFVRS IPDGSEKGFDF IALDLGGSSF RILRVQVNHE KNQNVHMESE VYDTPENIVH GSGSQLFDHV AECLGDFMEK RKIKDKKLPV GFTFSFPCQQ SKIDEAILIT WTKRFKASGV EGADVVKLLN KAIKKRGDYD ANIVAVNDT VGTMMTCGYD DQHCEVGLII GTGTNACYME ELRHIDLVEG DEGRMCINTE WGAFGDDGSL EDIRTEFDRE IDRGLNPGK QLFKEMVSGM YLGELVRLIL VKMAKEGLLF EGRITPELLT RGKFNSTSDVS AIEKNKEGLH NAKEILTRLG VEPSDDDCVS VQHVCTIVSF RSANLVAATL GAILNRLRDN KGTPRLRTTV GVDGSLYKTH PQYSRRFHKT LRRLVPDSDV RFLLESSESGS KGAAMVTAVA YRLAEQHRQI EETLAHFHILT KDMILLEVKKR MRAEMELGLR KQTHHNAVVK MLPSFVRRTP DGTENGDFLA LDLGGTNRV LLVKIRSGKK RTVEMHNKIY AIPIEIMQGT GEELFDHIVS CISDFLDYMG IKGPRMPLGF TFSFPCQQTS LDAGILITWT KGFKATDCVG HDVVTLRDA IKRREEFDLD VVAVVNDTVG TMMTCAYEEP TCEVGLIVGT GSNACYMEEM KNVEMVEGDQ GQMCINMEWG AFGDNGCLDD IRTHYDRLVD EYSLNAGKQR YEKMISGMYL GEIVRNILID FTKKGFLFRG QISETLKTRG IFETKFLSQI ESDRLALLQV RAILQQLGLN STCDDSIKLVK TVCGVVSRRR AQLCGAGMAA VVDKIRENRG LDRLNVTVG V DGTLYKLHPH FSRIMHQTVK ELSPKCNVSF LLEDGSGGK AALITAVGVR LRTEASS

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human Hexokinase 1 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.