

Human POLR2C Recombinant Protein

Catalog Number: PROTP19387

Overview

Product Name	Human POLR2C Recombinant Protein
Description	Human POLR2C Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-275aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 10% glycerol, 1mM DTT, 0.1M NaCl

Concentration

0.5mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, P19387, 1-275aa; MPYANQP TVRITELTDE NVKFIIENTD LAVANSIRRV FIAEVPIIAI DWVQIDANSS VLHDEFIAHR LGLIPLISDD
IVDKLQYSRD CTCEEFCPEC SVEFTLDVRC NEDQTRHVTS RDLISNSPRV IPVTSRNRDN DPNDYVEQDD ILIVKLRKGQ ELRLRAYAKK
GFGKEHAKWN PTAGVAFEYD PDNALRHTVY PKPEEWPKE YSELDEDESQ APYDPNGKPE RFYYNVESCGLRPETIVLS ALSGLKKKLS
DLQTQLSHEI QSDVLTIN

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human POLR2C Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.