

Human sTNF RI/TNFRSF1A Recombinant Protein

Catalog Number: PROTP19438-5

Overview

Product Name	Human sTNF RI/TNFRSF1A Recombinant Protein
Description	Recombinant human sTNF RI/TNFRSF1A protein expressed in E.coli with His-Tag. Sequence domain: 22-211aa. Suitable applications: SDS-PAGE, Denatured. Endotoxin not measured.
Source	E.coli
Physical Appearance	Liquid
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 10% glycerol

Storage

Can be stored at +2C to +8C for 1 week. For long term storage, aliquot and store at -20C to -80C. Avoid repeated freezing and thawing cycles.

Purity

> 85% by SDS-PAGE

USAGE

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.

