

Human recombinant ABP1/AOC1 protein

Catalog Number: PROTP19801-1

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Human recombinant ABP1/AOC1 protein
Description	Human recombinant ABP1/AOC1 protein expressed in HEK293 with His-Tag. Sequence domain: 20-751aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	HEK293
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

0.25mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, P19801, 20-751aa; EPSPGTLPRK AGVFSDLSNQ ELKAVHSFLW SKKELRLQPS STTTMAKNTV FLIEMLLPKK YHVLRFDDKG
ERHPVREARA VIFFGDQEHF NVTEFAVGPL PGPCYMRALS PRPGYQSSWA SRPISTAEYA LLYHTLQEAT KPLHQFFLNT TGFSFQDCHD
RCLAFTDVAP RGVASGQRRS WLIIQRYVEG YFLHPTGLEL LVDHGSTDAG HWAVEQVWYN GKFGSPPEEL ARKYADGEVD VVLEDPLPG
GKGHDSTEEL PLFSSHKPRG DFPSPHVSF PRLVQPHGPR FRLEGNAVLY GGWSFAFRLR SSSGLQVLNV HFGGERIAYE VSVQEAVLY

GGHTPAGMQT KYLDVGWGLG SVTHELAPGI DCPETATFLD TFHYDADDP VHYPRALCLF EMPTGVPLRR HFNSNFKGGF NFYAGLKGQV
LVLRTTSTVY NYDYIWDFIF YPNGVMEAKM HATGYVHATF YTPGLRHGT RLHTHLIGNI HTHLVHYRVD LDVAGTKNSF QTLQMKLENI
TNPWSPRHRV VQPTLEQTQY SWERQAAFRF KRKLPKYLLF TSPQENPWGH KRTYRLQIHS MADQVLPPGW QEEQAITWAR YPLAVTKYRE
SELCSSSIYH QNDPWHPV V FEQFLHNNEN IENEDLVAWV TVGFLHIPHS EDIPNTATPG NSVGFLLRPF NFFPEDPSLA SRDTVIVWPR
DNGPNYVQRW IPEDRDCSMP PPFSYNGTYR PV

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

[Human recombinant ABP1/AOC1 protein](#)

For Research Use Only. Not for use in diagnostic procedures.