

Mouse Hspa5 Recombinant Protein

Catalog Number: PROTP20029

Overview

Product Name	Mouse Hspa5 Recombinant Protein
Description	Mouse Hspa5 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 20-655aa. Application(s): SDS-PAGE.
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

0.5mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Mouse, P20029, 20-655aa; EEEDKKE DVGTVVIGIDL GTTYSCVGVF KNGRVEIIAN DQGNRITPSY VAFTPEGERL IGDAAKNQLT SNPENTVFDA KRLIGRTWND PSVQQDIKFL PFKVVEKTK PYIQVDIGGG QTKTFAPEEI SAMVLTKMKE TAEAYLGKKV THAVVTPAY FNDAQRQATK DAGTIAGLNV MRIINEPTAA AIAYGLDKRE GEKNILVFDL GGGTFDVSLL TIDNGVFEVV ATNGDTHLGG EDFDQRVMEH FIKLYKKKTG KDVRKDNRAV QKLRREVEKA KRALSSQHQA RIEIESFFEG EDFSETLTRA KFEELNMDLF RSTMKPVQKV LEDSDLKKS IDEIVLVGGS TRIPKIQQLV KEFFNGKEPS RGINPDEAVA YGAAVQAGVL SGDQDTGDLV LLDVCPLTLG IETVGGVMTK LIPRNTVVPT KKSQIFSTAS DNQPTVTIKV YEGERPLTKD NHLLGTFDLT GIPPAPRGVP QIEVTFEIDV NGILRVTAED KGTGNKNKIT ITNDQNRLTP EEIERMVNDA EKFAEEDKKL KERIDTRNEL ESYAYSLKNQ IGDKEKLGKG LSEDKETME KAVEEKIEWL ESHQDADIED FKAKKKELEE IVQPIISKLY GSGGPPPTGE EDTSEKDEL

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Mouse Hspa5 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.