

Human recombinant IBSP/Sialoprotein II protein

Catalog Number: PROTP21815-2

Overview

Product Name	Human recombinant IBSP/Sialoprotein II protein
Description	Human recombinant IBSP/Sialoprotein II protein expressed in HEK293 with His-Tag. Sequence domain: 17-317aa. Application(s): SDS-PAGE, Bioactivity. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	HEK293
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

0.5mg/ml (determined by Absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, P21815, 17-317aa; FSMKNLHRRV KIEDSEENG VFKYRPRYYLY KHAYFYPHLK RFPVQGSSDS SEENGDDSS EEEEEETS N
EGENNEESNE DEDSEAENTT LSATTLGYGE DATPGTGYTG LAAIQLPKKA GDITNKATKE KESDEEEEE EEGNENESE AEVDENEQGI
NGTSTNSTEA ENGNSSGGD NGEEGEEESV TGANAEDTTE TGRQKGTSK TTTSPNGGFE PTTTPQVYRT TSPFPGKTTT VEYEGEYET
GANEDNGYE IYSENGEPR GDNRYAYEDE YSYFKGGYD GYDGGNYYHH Q

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human recombinant IBSP/Sialoprotein II protein

For Research Use Only. Not for use in diagnostic procedures.