

Human PCMT Recombinant Protein

Catalog Number: PROTP22061

Overview

Product Name	Human PCMT Recombinant Protein
Description	Human PCMT Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-227aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 10% glycerol 0.1M NaCl

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Human, P22061, 1-227aa; MAWK SGGASHSELI HNLRKNGIIK TDKVFEVMLA TDRSHYAKCN PYMDSPQSIG FQATISAPHM HAYALELLFD
QLHEGAKALD VGSGSGILTA CFARMVGCTG KVIDIGDHIKE LVDDSINNVR KDDPTLLSSG RVQLVVG DGR MGYAEEAPYD AIHVGAAPV
VPQALIDQLK PGGRLILPVG PAGGNQMLEQ YDKLQDGSIK MKPLMGVIYV PLTDKEKQWS RWK

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human PCMT Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.