

## Human CCL1/I-309 Recombinant Protein

Catalog Number: PROTP22362

### Overview

Product Name	Human CCL1/I-309 Recombinant Protein
Description	Human CCL1/I-309 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 24-96aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 7.5) containing 2mM DTT, 10% glycerol, 50mM NaCl

### Concentration

0.5mg/ml (determined by Bradford assay)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles.

### Purity

> 85% by SDS-PAGE

### Amino Acid Sequence

Human, P22362, 24-96aa; MGSSHHHHHH SSGLVPRGSH MKSMQVPFSR CCFSFAEQEI PLRAILCYRN TSSICSNEGL IFKLKRGKEA  
CALDTVGWVQ RHRKMLRHCP SKRK [?/description\\_after\\_attributes?](#)

## Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human CCL1/I-309 Recombinant Protein