

Human TUBG1 Recombinant Protein

Catalog Number: PROTP23258

Overview

Product Name	Human TUBG1 Recombinant Protein
Description	Human TUBG1 Recombinant Protein expressed in Baculovirus with His-tag. Sequence domain: 1-451aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	Baculovirus
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 40% glycerol, 0.1M NaCl, 2mM DTT, 50mM imidazole.

Concentration

0.25mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 85% by SDS-PAGE

Amino Acid Sequence

Human, P23258, 1-451aa; MPREIITLQL GQCGNQIGFE FWKQLCAEHG ISPEGIVEEF ATEGTRKDV FFYQADDEHY IPRAVLLDLE PRVIHSILNS
PYAKLYNPEN IYLSEHGGGA GNNWASGFSQ GEKIHEDIFD IIDREADGSD SLEGFVLCBS IAGGTGSLG SYLLERLNDR YPKKLVQTYS
VFPNQDEMSD VVVQPYNLL TLKRLTQNAD CVVLDNTAL NRIATDRLHI QNPSFSQINQ LVSTIMSAST TTLRYPGYMN NDLIGLIASL
IPTPRLHFLM TGYTPLTTDQ SVASVRKTTV LDVMRLLQP KNVMVSTGRD RQTNHCYIAI LNIIQGEVDP TQVHKSLQRI RERKLANFIP
WGPASIQVAL SRKSPYLPSSA HRVSGLLMMAN HTSISSLFER TCRQYDKLRK REAFLEQFRK EDMFKDNFDE MDTSREIVQQ LIDEYHAATR
PDYISWGTQE Q

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human TUBG1 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.