

Human WARS Recombinant Protein

Catalog Number: PROTP23381

Overview

Product Name	Human WARS Recombinant Protein
Description	Human WARS Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-471aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 0.1M NaCl, 1mM DTT, 10% glycerol

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, P23381, 1-471aa; MGSSHHHHHH SGLVPRGSH MPNSEPASLL ELFNSIATQG ELVRSKAGN ASKDEIDSAV KMLVSLKMSY
KAAAGEDYKA DCPGPNPAPT SNHGPDTEA EEDFVDPWTV QTSSAKGIDY DKLIVRFSS KIDKELINRI ERATGQRPHH FLRRGIFFSH
RDMNQVLDAY ENKKPFYLYT GRGPSSEAMH VGHLIPFIFT KWLQDVFNVP LVIQMTDDEK YLWKDLTLDQ AYSYAVENAK DIIACGFDIN
KTFIFSDLDY MGMSSGFYKN VVKIQKHVTF NQVKGIFGFT DSDCIGKISF PAIQAAAPSFS NSFPQIFRDR TDIQCLIPCA IDQDPYFRMT
RDVAPRIGYP KPALLHSTFF PALQGAQTKM SASDPNSSIF LTDTAKQIKT KVNKHAFSSG RDTIEEHRQF GGNCDVDVSF MYLTFLEDD
DKLEQIRKDY TSGAMLTGEL KKALIEVLQP LIAEQARRK EVTDEIVKEF MTPRKLSFDF Q

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human WARS Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.