

Mouse Th Recombinant Protein

Catalog Number: PROTP24529

Overview

Product Name	Mouse Th Recombinant Protein
Description	Mouse Th Recombinant Protein expressed in Baculovirus with His-tag. Sequence domain: 1-498aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

0.25mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Mouse, P24529, 1-498aa; MPTPSAS SPQPKGFERR VSEQDTKQAE AVTSPRFIGR RQSLIEDARK EREAAAAAAAA AAVASAEPGN PLEAVVFEER
DGNVNLNLLF SLRGTKPSSL SRALKVFETF EAKIHHLLETR PAQRPLAGSP HLEYFVRFEV PSGDLAALLS SVRRVSDDV R SAREDKVPWF
PRKVSELDKC HHLVTKFDPD LDLDHPGFSD QAYRQRRKLI AEIAFQYKQG EPIPHVEYTK EEIATWKEVY ATLKGLYATH ACREHLEAFQ
LLERYCGYRE DSIPQLEDVS HFLKERTGFQ LRPVAGLLSA RDFLASLAFR VFQCTQYIRH ASSPMHSPEP DCHELLGHV PMLADRTFAQ
FSQDIGLASL GASDEEIEKL STVYWFTVEF GLCKQNGELK AYGAGLLSSY GELLHSLSEE PEVRAFD PDT AAVQPYQDQT YQPVYFVSES
FSDAKDKLRN YASRIQRPFS VKFDPYTLAI DVLDSPHTIR RSLEGVQDEL HTLTQALSAI S

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Mouse Th Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.