

## Rat ALT1/GPT1 Recombinant Protein

Catalog Number: PROTP25409

### Overview

Product Name	Rat ALT1/GPT1 Recombinant Protein
Description	Rat ALT1/GPT1 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-496aa. Application(s): Enzyme Activity, SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

### Concentration

1mg/ml (determined by absorbance at 280nm)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 90% by SDS-PAGE

### Amino Acid Sequence

Rat, P25409, 1-496aa; MGSSHHHHHH SGLVPRGSH MGSMASRVND QSQASRNLK GKVLTLDTMN PCVRRVEYAV RGPVQRAL  
LEQLRQGVK KPFTVIRAN IGDAQAMGQR PITFFRQVLA LCVYPNLLSS PDFPEDAKRR AERILQACGG HSLGAYSISS GIQPIREDVA  
QYIERRDGGI PADPNNIFLS TGASDAIVTM LKLLVSGEGR ARTGVLIPIP QYPLYSALA ELDAVQVDYY LDEERAWALD IAE LRALCQ  
ARDRCCPRVL CVINPGNPTG QVQTRECIEA VIRFAFKEGL FLMADEVYQD NVYAEGSQFH SFKKVLMEMG PPySTQQELA SFHSVSKGYM  
GECGFRGGYV EVVNMDAEVQ KQMGKLMsvr LCPPVPGQAL MDMVVSPPTP SEPSFKQFQA ERQEVLAELA AKAKLTEQVF NEAPGIRCNP  
VQGAMYSFPQ VQLPLKAVQR AQELGLAPDM FFCLCLEET GICVVPVSGF GQQEGTYHFR MTILPPMEKL RLLLEKLSHF HAKFTHEYS

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Rat ALT1/GPT1 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.