

Human BRD2 Recombinant Protein

Catalog Number: PROTP25440

Overview

Product Name	Human BRD2 Recombinant Protein
Description	Human BRD2 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-455aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	Phosphate-Buffered Saline (pH 7.4) containing, 10% glycerol, 1mM DTT

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, P25440, 1-455aa; MLQNVTP HNKLPGEGNA GLLGLGPEAA APGKRIRKPS LLYEGFESPT MASVPALQLT PANPPPPEVS NPKKPGRVTN QLQYLHKVVM KALWKHQFAW PFRQPVDVAVK LGLPDYHKII KQPMDMGTIK RRLENNYYWA ASECMQDFNT MFTNCYIYNK PTDDIVLMAQ TLEKIFLQKV ASMPQEEQEL VVTIPKNSHK KGAKLAALQG SVTSAHQVPA VSSVSHTALY TPPPEIPTTV LNIPHPSVIS SPLKSLHSA GPPLAVTAA PPAQPLAKKK GVKRKADTTT PTPTAILAPG SPASPPGSLE PKAARLPPMR RESGRPIKPP RKDLPDSQQQ HQSSKKGKLS EQLKHCNGIL KELLSSKHAA YAWPFYKPVD ASALGLHDYH DIIKHPMDLS TVKRKMNENR YRDAQEFAAD VRLMFSNCYK YNPPDHDVVA MARKLQDVFE FRYAKMPD

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human BRD2 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.