

## Human I $\kappa$ B-alpha Recombinant Protein

Catalog Number: PROTP25963

### Overview

Product Name	Human I $\kappa$ B-alpha Recombinant Protein
Description	Human I $\kappa$ B-alpha Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-317aa. Application(s): SDS-PAGE.
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 1mM DTT, 20% glycerol, 0.1M NaCl

### Concentration

0.5mg/ml (determined by Bradford assay)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 85% by SDS-PAGE

### Amino Acid Sequence

Human, P25963, 1-317aa; MGSSHHHHHH SGLVPRGSH MFQAAERPQE WAMEGPRDGL KKERLLDDRH DSGLD SMKDE EYEQMVKELQ  
EIRLEPQEVV RGSEPWKQQL TEDGDSFLHL AIIHEEKALT MEVIRQVKGD LAFLNFQNNL QQTPLHLAVI TNQPEIAEAL LGAGCDPELR  
DFRGNTPLHL ACEQGCLASV GVL TQSC TTP HLH SILKATN YNGHTCLHLA SIHG YLGIVE LLVSLGADVN AQEPCNGRTA LHLAVDLQNP  
DLVSLLLKCG ADVNRVTYQG YSPYQLTWGR PSTRIQQQLG QLTLENLQML PESEDEESYD TESEFTEFTE DELPYDDCVF GGQRLTL

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Human I $\kappa$ B-alpha Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.