

Human CNTF Recombinant Protein

Catalog Number: PROTP26441

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Human CNTF Recombinant Protein
Description	Human CNTF Recombinant Protein expressed in E. coli (non-tagged). Sequence domain: 1-200aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	Non-Tagged
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 1mM DTT, 10% glycerol

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, P26441, 1-200aa; MAFTEHSPLT PHRRDLCRS IWLARKIRSD LTALTESYVK HQGLNKNINL DSADGMPVAS TDQWSELTEA
ERLQENLQAY RTFHVLLARL LEDQQVHFTP TEGDFHQAIH TLLQVAAFA YQIEELMILL EYKIPRNEAD GMPINVGDDG LFEKKLWGLK
VLQELSQWTV RSIHDLRFIS SHQTGIPARG SHYIANNKKM

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Human CNTF Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.