

Human HMGB2 Recombinant Protein

Catalog Number: PROTP26583

Overview

Product Name	Human HMGB2 Recombinant Protein
Description	Human HMGB2 Recombinant Protein expressed in Baculovirus with His-tag. Sequence domain: 1-209aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	Baculovirus
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 1mM DTT, 30% glycerol, 0.1M NaCl

Concentration

0.5mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, P26583, 1-209aa; MGKGDPN KPRGKMSSYA FFVQTCREEH KKKHPDSSVN FAEFSKKCSE RWKTMSAKEK SKFEDMAKSD
KARYDREMKN YVPPKGDKKG KKKDPNAPKR PPSAFFLFCSE EHRPKIKSEH PGLSIGDTAK KLGEMWSEQS AKDKQPYEQK AAKLKEKYEK
DIAAYRAK GK SEAGKKG PGR PTGSKKKNEP EDEEEEEEE DEDEEEDEE EE

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human HMGB2 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.