

Mouse recombinant IL-3R alpha/IL3RA protein

Catalog Number: PROTP26952

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Mouse recombinant IL-3R alpha/IL3RA protein
Description	Mouse recombinant IL-3R alpha/IL3RA protein expressed in Baculovirus with hlgG-His-Tag. Sequence domain: 17-331aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	hlgG-His-Tag
Form	Liquid
Source	Baculovirus
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 0.1M NaCl, 30% glycerol

Concentration

0.25mg/ml (determined by Absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Mouse, P26952, 17-331aa; SDLAAVREAP PTAVTTPIQN LHIDPAHYTL SWDPAPGADI TTGAFCRKGR DIFVWADPGL ARCSFQSLSL
CHVTNFTVFL GKDRVAGSI QFPDGGDH EAAQDLRCW VHEGQLSCQW ERGPKATGDV HYRMFWRDVR LGPAHNRECP HYHSLDVNTA
GPAPHGGHEG CTLDLDTVLG STPNPDLVP QVTITVNGSG RAGPVPCMDN TVDLQRAEVL APPTLTVECN GSEAHARWVA RNRFHGLLG
YTLQVNQSSR SEPQEYNVSI PHFWVPNAGA ISFRVKSRSE VYPRKLSSWS EAWGLVCPPE VMPVK

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

[Mouse recombinant IL-3R alpha/IL3RA protein](#)

For Research Use Only. Not for use in diagnostic procedures.