

## Mouse CD40 Ligand/TNFSF5 Recombinant Protein

Catalog Number: PROTP27548

### Overview

Product Name	Mouse CD40 Ligand/TNFSF5 Recombinant Protein
Description	Mouse CD40 Ligand/TNFSF5 Recombinant Protein expressed in Baculovirus with His-tag. Sequence domain: 112-260aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

### Concentration

0.25mg/ml (determined by absorbance at 280nm)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 95% by SDS-PAGE

### Amino Acid Sequence

Mouse, P27548, 112-260aa; MQRGDED PQIAAHVVSE ANSNAASVLQ WAKKGYTMMK SNLVMLENGK QLTVKREGLY YVYTQVTFCS NREPSSQRPF IVGLWLKPSS GSERILLKAA NTHSSSQLCE QQSVHLGGVF ELQAGASV FV NVTEASQVIH RVGFSSFGLL KL

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Mouse CD40 Ligand/TNFSF5 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.