

Human Calnexin Recombinant Protein

Catalog Number: PROTP27824

Overview

Product Name	Human Calnexin Recombinant Protein
Description	Human Calnexin Recombinant Protein expressed in E. coli (non-tagged). Sequence domain: 21-481aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	Non-Tagged
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 1mM DTT and 20% glycerol

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, P27824, 21-481aa; MHDGHDDDDVI DIEDDLDDVI EEVEDSKPDT TAPPSSPKVT YKAPVPTGEV YFADSFDRGT LSGWILSKAK
KDDTDDEIAK YDGKWEVEEM KESKLPDGDKG LVLMsRAKHH AISAKLNKPF LFDTKPLIVQ YEVENFQNGIE CGGAYVKLLS KPELNLDQF
HDKTPYTIMF GPKKCGEDYK LHFIFRHKNP KTGIIYEEKHA KRPDADLKTY FTDKKTHLYT LILNPDNSFE ILVDQSVVNS GNLLNDMTPP
VNPSREIEDP EDRKPEDWDE RPKIPDPEAV KPDDWDEDAP AKIPDEEATK PEGWLDDPEE YVPDPAEKP EDWDEDMDGE WEAPQIANPR
CESAPGCGVW QRPVIDNPNY KGKWKPPMID NPSYQGIWKP RKIPNPDFFE DLEPFRMTPF SAIGLELWSM TSDIFFDNFI ICADRRIVDD
WANDGWGLKK AADGAAEPGV VGQMIEAAEE RP

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human Calnexin Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.