

## Mouse Biglycan Recombinant Protein

Catalog Number: PROTP28653

### Overview

Product Name	Mouse Biglycan Recombinant Protein
Description	Mouse Biglycan Recombinant Protein expressed in Baculovirus with hlgG-His-tag. Sequence domain: 38-369aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	hlgG-His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

### Concentration

0.5mg/ml (determined by absorbance at 280nm)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 90% by SDS-PAGE

### Amino Acid Sequence

Mouse, P28653, 38-369aa; ADPDEEASGS DTTSGVPDL D SVTPTFSAMC PFGCHCHLRV VQCS DLGLKT VPKEISPD TT LLDLQNNDIS  
ELRKDDFKGL QHLYALVLVN NKISKIHEKA FSPLRKLQKL YISK NHLVEI PPNLPSSLVE LRIHDNRIRK VPKG VFSGLR NMNCIEMGGN  
PLENSGFEPG AFDGLKLNLYL RISEAKLTGI PKDLPETLNE LHL DHNKIQA IELEDLLRYS KLYRLGLGHN QIRMIENGSL SFLPTLRELH LDNNKLSRVP  
AGLPDLKLLQ VVYLHSNNIT KVGINDFCPM GFGVKRAYYN GISLFNNPVP YWEVQPATFR CVTDRLAIQF GNYKKLEPKS CDKTHTCPPC  
PAPELLGGPS VFLFPPKPKD TLMISRTPEV TCVVVDVSHE DPEVKFNWYV DGVEVHNAKT KPREEQYNST YRVVSVLTVL HQDWLNGKEY  
KCKVSNKALP APIEKTISKA KGQPREPQVY TLPPSRDELT KNQVSLTCLV KGFYPSDIAV EWESNGQPEN NYKTTTPVLD SDGSFFLYSK  
LTVDKSRWQQ GNVFSCSVMH EALHNHYTQK SLSLSPGKHH HHHH

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Mouse Biglycan Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.