

PGRN Progranulin Human Recombinant Protein

Catalog Number: PROTP28799

Overview

Product Name	PGRN Progranulin Human Recombinant Protein
Description	Progranulin Human Recombinant produced in HEK is a single, glycosylated, polypeptide chain containing 1-593 amino acids and having a molecular mass of 74kDa. The Progranulin is purified by standard chromatographic techniques.
Size	2ug, 10ug, 100ug
Tag	
Form	Sterile Filtered White lyophilized (freeze-dried) powder.
Source	HEK293
Formulation	The protein contains 1xPBS.

Concentration

Storage

Lyophilized Progranulin although stable at room temperature for 3 weeks, should be stored desiccated below -18°C. Upon reconstitution PGRN should be stored at 4°C between 2-7 days and for future use below -18°C. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Please prevent freeze-thaw cycles.

Purity

Greater than 98.0% as determined by SDS-PAGE.

!description_after_attributes! 2/description_after_attributes2

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



PGRN Progranulin Human Recombinant Protein