

Human CD38 Recombinant Protein

Catalog Number: PROTP28907

Overview

Product Name	Human CD38 Recombinant Protein
Description	Human CD38 Recombinant Protein expressed in Baculovirus with His-tag. Sequence domain: 43-300aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	Baculovirus
Formulation	50mM MES buffer (pH 5.0) containing 100mM NaCl, 10% glycerol

Concentration

1mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, P28907, 43-300aa; ADPEFVPRWR QQWSGPGTTK RFPETVLARC VKYTEIHPEM RHVDCQSVWD AFKGAFISKH PCNITEEDYQ
PLMKLGTQTV PCNKILLWSR IKDLAHQFTQ VQRDMFTLED TLLGYLADDL TWCGEFNTSK INYQSCP DWR KDCSNNPVSV FWKTVSRRFA
EAACDVVHVM LNGRSKIFD KNSTFGSVEV HNLQPEKVQT LEAWVIHGGR EDSRDL CQDP TIKELESII S KRNIQFSCKN IYRPDKFLQC
VKNPEDSSCT SEIH HHHHHH

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human CD38 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.