

## Rat recombinant Erythropoietin/EPO protein

Catalog Number: PROTP29676

### Overview

Product Name	Rat recombinant Erythropoietin/EPO protein
Description	Rat recombinant Erythropoietin/EPO protein expressed in HEK293 with His-Tag. Sequence domain: 27-192aa. Application(s): SDS-PAGE, Bioactivity. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	HEK293
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

### Concentration

1mg/ml (determined by Absorbance at 280nm)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 95% by SDS-PAGE

### Amino Acid Sequence

Rat, P29676, 27-192aa; APPRLIC DSRVLERYIL EAKEAENVTM GCAEGPRLSE NITVPDTKVN FYAWKRMKVE EQAVEVWQGL SLLSEAILQA QALQANSSQP PESLQLHIDK AISGLRSLTS LLRVLGAQKE LMSPPDATQA APLRRTLADT FCKLFRVYSN FLRGKCLKLYT GEACRRGDR

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Rat recombinant Erythropoietin/EPO protein

For Research Use Only. Not for use in diagnostic procedures.