

Human recombinant CD6 protein

Catalog Number: PROTP30203

Overview

Product Name	Human recombinant CD6 protein
Description	Human recombinant CD6 protein expressed in Baculovirus with hIgG-His-Tag. Sequence domain: 18-402aa. Application(s): SDS-PAGE, Bioactivity. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	hIgG-His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

0.5mg/ml (determined by Absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

90% by SDS-PAGE

Amino Acid Sequence

Human, P30203, 18-402aa; HPSPAPP DQLNTSSAES ELWEPGERLP VRLTNGSSSC SGTVEVRLEA SWEPACGALW DSRAAEAVCR ALGCGGAEAA SQLAPPTPEL PPPPAAGNTS VAANATLAGA PALLCSGAEW RLCEVVEHAC RSDGRRARVT CAENRALRLV DGGGACAGRV EMLEHGEWGS VCDDTWDLED AHVVCRLGCG GWAVQALPGL HFTPGRGPIH RDQVNCSGAE AYLWDCPGLP GQHYCGHKED AGAVCSEHQSWRLTGGADRC EGQVEVHFRG VVNTVCDSEW YPSEAKVLCQ SLGCGTAVER PKGLPHSLSG RMYYSNCGEE LTLNCSWRF NNSNLCSQSL AARVLCASR SLHNLSTPEV PASVQVTIE SSVTVKIENK ESREMLLL

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human recombinant CD6 protein

For Research Use Only. Not for use in diagnostic procedures.