

## Human CDC25A Recombinant Protein

Catalog Number: PROTP30304

### Overview

Product Name	Human CDC25A Recombinant Protein
Description	Human CDC25A Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-524aa. Application(s): SDS-PAGE.
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 5mM DTT, 20% glycerol, 0.2M NaCl, 1mM EDTA

### Concentration

0.25mg/ml (determined by Bradford assay)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 80% by SDS-PAGE

### Amino Acid Sequence

Human, P30304, 1-524aa; MRGSHHHHHH GMASMTGGQQ MGRDLYDDDD KDRWGSMELG PEPPHRRRL FACSPPPASQ PVVKALFGAS  
AAGGLSPVTN LVTMDQLQG LGSDYEQPLE VKNNSNLQRM GSSESTDSGF CLDSPGPLDS KENLENPMRR IHSLPQKLLG CSPALKRSHS  
DSLHDHIFQL IDPDENKENE AFEFKKVPVR VSRGCLHSHG LQEGKDLFTQ RQNSAPARML SSNERDSSEP GNFIPLFTPQ SPVTATLSDE  
DDGFVDLLDG ENLKNEETP SCMASLWTAP LVMRTTNLDN RCKLFDSPSL CSSSTRSVLK RPERSQEES PGSTKRRKSM SGASPKESTN  
PEKAHETLHQ SLSLASSPKG TIENILDNDP RDLIGDFSKG YLFHTVAGKH QDLKYISPEI MASVLNGKFA NLIKEFVIID CRYPYEYEGG  
HIKGAVNLHM EEEVEDFLLK KPIVPTDGKR VIVVFHCFES SERGPRMCRY VRERDRLGNE YPKLHYPELY VLKGGYKEFF MKCQSYCEPP  
SYRPMHHEDF KEDLKKFR TK SRTWAGEKSK REMYSRLKKL

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human CDC25A Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.