

Human Cyclophilin D/Cyclophilin F/PPIF Recombinant Protein

Catalog Number: PROTP30405

Overview

Product Name	Human Cyclophilin D/Cyclophilin F/PPIF Recombinant Protein
Description	Human Cyclophilin D/Cyclophilin F/PPIF Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 30-207aa. Application(s): SDS-PAGE, Enzyme Activity. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 7.5) containing 1mM DTT, 10% glycerol

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Human, P30405, 30-207aa; MGSSHHHHHH SGLVPRGSH MCSKSGSDPS SSSSSGNPLV YLDVDANGKP LGRVVLELKA DVVPKTAENF RALCTGEKGF GYKGFSTFHRV IPSFMCQAGD FTNHNGTGGK SIYGSRFPDE NFTLKHVGPV VLSMANAGPN TNGSQFFICT IKTDWLDGKH VVFGHVKEGM DVVKKIESFG SKSGRTSKKI VITDCGQLS

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human Cyclophilin D/Cyclophilin F/PPIF Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.