

## Human NMT1 Recombinant Protein

Catalog Number: PROTP30419

### Overview

Product Name	Human NMT1 Recombinant Protein
Description	Human NMT1 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-496aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	Phosphate-Buffered Saline (pH 7.4)

### Concentration

1mg/ml (determined by Bradford assay)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 85% by SDS-PAGE

### Amino Acid Sequence

Human, P30419, 1-496aa; MADESET AVKPPAPPLP QMMEGNGNGH EHCSDCENEE DNSYNRGGLS PANDTGAKKK KKKQKKKKEK  
GSETDSAQDQ PVKMNSLP AE RIQEIQKAI E LFSVGQGPAK TMEEASKRSY QFWDTPVVPK LGEVVNTHGP VEPDKDNIRQ EPYTL PQGFT  
WDALDLGDRG VLKELYLLN ENYVEDDDNM FRFDYSPEFL LWALRPPGWL PQWHCGVRVV SSRKLVGFIS AIPANIHIYD TEKKMVEINF  
LCVHKKLRSK RVAPVLIREI TRRVHLEGIF QAVYTAGVVL PKPVGTCRYW HRSLNPRKLI EVKFSHLSRN MTMQR TMKLY RLPETPKTAG  
LRPMETKDIP VVHQLTRYL KQFHLTPVMS QEEVEHWFYP QENIIDTFVV ENANGEVTDF LSFYTL PSTI MNHPH THKSLK AAYSFYNVHT  
QTPLLDLMSD ALVLAKMKGF DVFNALDLME NKTFLKLF GIGDGNLQYY LYNWKCP SMG AEKVGLVLQ

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Human NMT1 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.