

Rantes Mouse Recombinant Protein (CCL5)

Catalog Number: PROTP30882

Overview

Product Name	Rantes Mouse Recombinant Protein (CCL5)
Description	Rantes Mouse Recombinant produced in E. coli is a single, non-glycosylated polypeptide chain containing 68 amino acids and having a molecular mass of 7876 Dalton. The Mouse Rantes is purified by proprietary chromatographic techniques.
Size	5ug, 20ug, 1mg
Tag	
Form	Sterile Filtered White lyophilized (freeze-dried) powder.
Source	Escherichia coli
Formulation	Lyophilized from a concentrated (1mg/ml) solution in water containing no additives.

Concentration

Storage

Lyophilized Mouse Rantes although stable at room temperature for 3 weeks, should be stored desiccated below -18°C. Upon reconstitution Mouse CCL5 should be stored at 4°C between 2-7 days and for future use below -18°C. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Please prevent freeze-thaw cycles.

Purity

Greater than 95.0% as determined by (a) Analysis by RP-HPLC and (b) Analysis by SDS-PAGE.

Amino Acid Sequence

SPYGSDTTPC CFAYLSLALP RAHVKEYFYT SSKCSNLAVV FVTRRRNRQVC ANPEKKWVQE YINYLEMS

[?/description_after_attributes?](#)

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Rantes Mouse Recombinant Protein (CCL5)